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## General Warnings and Safety Instructions

### **You Have Chosen a Hearing System from HANSATON.**

We take pleasure in your choice and congratulate you on your decision, because with a HANSATON hearing system you own a high quality product of the most modern technology.

For 60 years, we have been helping people with our devices to achieve better hearing and better speech understanding, thus a higher quality of life.

Behind our products stand many years of development and production experience, as a corporate philosophy centered on people and their needs.

If you read this brochure carefully and follow the advice on use and care, you will get much enjoyment from your new hearing system.

This product has passed every phase of the independently-certified HANSATON quality management system.

### **Intended Use.**

Hearing systems are intended to compensate hearing loss of hearing-impaired persons.

Diagnosis of hearing loss, prescription and fitting of a hearing system must be performed by hearing health specialists (e.g., ENT specialists,

audiologists, or hearing system professionals). The consulting must also be performed by a hearing health specialist.

Read and follow the instructions of this user guide to avoid injuries or damage to the hearing system.

Hearing systems will not restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions. The use of hearing systems is only part of hearing rehabilitation; auditory training and lip reading instruction may be required as well.

In most cases, infrequent use of hearing systems does not provide full benefit. Once you have become accustomed to your hearing systems, wear your hearing systems everyday all day.

### Symbols.

With the CE symbol, Hansaton confirms that this Hansaton product – including accessories – meets the requirements of the Medical Devices Directive 93/42/EEC as well as the Radio Equipment Directive 2014/53/EU on radio and telecommunications transmitters. The numbers after the CE symbol are the numbers of the consulted certified bodies under the above mentioned directives.

This symbol indicates that it is important for the user to read and take into account the relevant information in this user guide.

This symbol indicates that it is important for the user to pay attention to the relevant warning notices in this user guide.

Important information for handling and effective use of the product.

This symbol indicates that the products described in these user instructions adhere to the requirements for an application part of Type B of EN 60601-1. The surface of the hearing system is specified as application part of Type B.

Copyright symbol

This symbol shall be accompanied by the name and the address of the authorised representative in the European Community.

This symbol shall be accompanied by the name and the address of the manufacturer (who are placing this device on the market).

The symbol with the crossed-out garbage bin is to make you aware that this device may not be thrown away as normal household waste. Please dispose of old or unused devices, at waste disposal sites intended for electronic waste, or give your device to your hearing healthcare professional for disposal. Proper disposal protects the environment and health.

 **Bluetooth®** The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Hansaton is under license. Other trademarks and trade names are those of their respective owners.

### Warnings.

- Do not use your hearing aid in potentially explosive environments.
- If you experience dizziness, drainage or bleeding from the ear, or sudden hearing loss, see your Physician immediately.
- A hearing aid may cause an accelerated accumulation of cerumen (earwax).
- Keep hearing aids out of the reach of children, pets, and individuals with diminished mental capacity.
- Hearing aids should be used only as directed and adjusted by your trained hearing care professional.
- A hearing aid may stop functioning, e.g. if the battery goes dead. You should be aware of this possibility, in particular if you are in traffic or otherwise depend on warning sounds.
- Do not switch hearing aids between ears or allow others to wear or use your hearing aids.

The following is only applicable for persons with active implantable medical devices (i.e. pacemakers, defibrillators, etc.):

- Keep the wireless hearing aid at least 15 cm (6 inches) away from the active implant.
- If you experience any interference, do not use the wireless hearing aids and contact the manufacturer of the active implant. Please note that interference can also be caused by power lines, electrostatic discharge, airport metal detectors etc.

## TINNITUS MANAGER

### Important Safety Information.

The Tinnitus Manager is a broadband sound generator. It provides a means of sound enrichment therapy that can be used as part of a personalized tinnitus management program to provide temporary relief from tinnitus.

The underlying principle of sound enrichment is to provide supplementary noise stimulation which can help defocus your attention from your tinnitus and avoid negative reactions. Sound enrichment, coupled with instructional counseling, is an established approach to managing tinnitus.

Your hearing aid is equipped with a Tinnitus Manager. The use of the Tinnitus Manager should be only on the advice and in consultation with your Hearing Care Professional.

### Warnings about the Tinnitus Manager.

 **WARNING!** Risk of further impairment to the user's hearing health. There are some potential concerns associated with the use of any sound generating tinnitus therapy device. Discontinue use and seek medical evaluation if any of the following conditions occur:

- Chronic skin irritation on, near, or around the site of device placement.
- Unusual side effects (e.g. dizziness, nausea, headaches, heart palpitations).
- Perceived decrease in auditory function (e.g. decreased loudness, speech not as clear).

 **WARNING!** Risk of further impairment to the user's hearing health. The volume of the Tinnitus Manager can be set to a level which could lead to permanent hearing damage when used for a prolonged period of time.

- Tinnitus Manager should never be used at uncomfortable levels.

### For Hearing Care Professionals:

The target population is primarily the adult population over 18 years of age. The patient may have some control of the level or volume of the signal and the patient should discuss this adjustment as well as his or her comfort level and sound of the signal with their Hearing Care Professional.

 **WARNING!** Risk of further impairment to the user's hearing health. The volume of the Tinnitus Manager can be set to a level which could lead to permanent hearing damage when used for a prolonged period of time.

- Should the Tinnitus Manager be set to such a level in the hearing aid, the maximum safe duration per day must be adapted to the sound level (dBA) as shown in the following table (29 CFR 1910.95 by the Occupational Safety and Health Administration (OSHA)).

In this case, the health care professional must inform the user of the device accordingly.

D

EN

Sound level dBA slow response	Maximum Safe Duration per day
82 dB (A)	24 hours
85 dB (A)	16 hours
90 dB (A)	8 hours

F

- The Tinnitus Manager should never be used at uncomfortable levels.

E

 **WARNING!** Risk of further impairment to the user's hearing health.

A Hearing Care Professional should advise a prospective tinnitus noiser user to consult promptly with a licensed physician (preferably an ear specialist) before using the Tinnitus Manager if the Hearing Care Professional determines through inquiry, actual observation, or review or any other available information concerning the prospective user that the prospective user has any of the following conditions:

- Visible congenital or traumatic deformity of the ear.
- History of active drainage from the ear within the previous 90 days.
- History of sudden or rapidly progressive hearing loss within the previous 90 days.
- Acute or chronic dizziness.
- Unilateral hearing loss of sudden or recent onset within the previous 90 days.

### General Notes on Safety.

**Your hearing system is sensitive to extreme heat, high humidity, strong magnetic fields (>0.1T), Xrays, and mechanical stress.**

- Do not expose your hearing system to extreme temperature or high humidity.
- Do not leave it in direct sunlight.
- Do not wear it in the shower or when you apply makeup, perfume, aftershave, hairspray, or suntan lotion.
- Do not wear your hearing system when you are exposed to short-waves, a strong magnetic field, a high frequency field, or X-rays.
- Do not place your hearing systems in the microwave.

When using mobile telephones, however, electromagnetic-field effects cannot be ruled out.

Your hearing systems operate in the frequency range of 2.4 GHz–2.48 GHz. When flying please check if flight operator requires devices to be switched into flight mode.

Allergic reactions to hearing systems are unlikely. However, if you experience itching, redness, soreness, inflammation or a burning sensation in or around your ears, inform your hearing healthcare professional and contact your physician.

In the unlikely case that any parts remain in the ear canal after the

removal of the hearing system, contact a physician immediately.

 **WARNING! Choking hazard posed by small parts.**

- Keep hearing systems, batteries, and accessories out of children's reach.
- If swallowed, immediately consult a physician or hospital.
- If infants, small children, or mentally disabled persons need to wear hearing systems, ensure adequate supervision.

 **WARNING! Risk in explosive environments (e.g., mining areas)!**

- Don't use your hearing systems (and remote control) in areas where there is a danger of explosions.

In some countries restrictions on the use of wireless equipment exist.

- Consult local authorities for further information.

 **WARNING!**

- Never use a tip that has been worn by someone else or if it becomes dirty, torn, or disfigured.
- Only use original tips. Tips from other manufacturers will not yield the best performance or proper fit.
- Do not insert the hearing system if you think you may have a foreign object in the ear canal.
- The hearing system and tips are made of modern non-allergenic materials. Still, in rare cases, skin irritation can occur. If you experience discomfort or skin irritation in the ear or in the ear canal, discontinue use, and contact your Physician immediately.
- If tip is lodged in your ear canal, please see a Physician immediately to ensure it is safely removed.

### Safety Notes for Batteries.

 **WARNING! Never attempt to recharge batteries that are not designed for that purpose (danger of explosion)!**

**Leaking batteries damage the system.**

- Remove the battery when the system is not in use for a prolonged period of time.
- To preserve the battery, turn your hearing systems off when not in use.

**Batteries contain harmful substances that pollute the environment.**

- Do not dispose of batteries in household trash.
- Dispose of batteries according to federal and local regulations, or return them to your hearing system professional.

## D Warnings about Magnets.

EN

 **WARNING! Magnets may interfere with the operation of medical devices (e.g., pacemakers and magnetic valves).**

F

- Use a magnet only if not in doubt about the safe distance.
- The safe distance between pacemaker and magnet should be at least 30 cm (12 inches).

E

 **WARNING!**

- Keep magnets away from children. If swallowed consult a physician or hospital immediately.

 **Magnets can disturb electrical devices and delete stored data.**

- Keep magnets away from computers, color monitors, television sets, CD, video and audio discs and other electronic equipment/devices.

When operating the hearing system in telecoil mode, the wireless function may cause a pulsing noise. This is due to the pick-up of the transmission signal.

- When your hearing system is in telecoil mode please use the remote control only at a distance greater than 30 cm (12 inches).

## Additional Warnings about Magnets Affixed to the Telephone.

- Be sure the magnet is securely affixed to the telephone.
- If the magnet falls into your ear, contact your hearing healthcare professional.
- Too high distortion during dialing or phoning may mean that the phone handset is stressed by the magnet. To avoid any damage, please move the magnet to another place on the telephone receiver.

## Domes and Receiver Molds.

- Always wear the Receiver Unit with the Dome.
- Make sure that the Dome is completely attached to the Receiver Unit.
- If the Dome remains in your ear, have it removed by a medical professional.
- Replace the Dome if it breaks.
- Insert the Receiver Unit with the Dome carefully and not too deeply into the ear.

## Note to Hearing Healthcare Professional:

Domes should never be fitted on patients with perforated eardrums, exposed middle ear cavities, or surgically altered ear canals. In the case of such a condition, we recommend to use a customized ear mold.

## Mobile phone compatibility rating.

Some hearing system users have reported a buzzing sound in their hearing systems when they are using mobile phones, indicating that the mobile phone and hearing system may not be compatible.

According to the ANSI C63.19 standard (ANSI C63.19-2011 American National Standard Methods of Measurement of Compatibility Between Wireless Communications Devices and Hearing Aids), the compatibility of a particular hearing aid and mobile phone can be predicted by adding the rating for the hearing aid immunity to the rating for the mobile phone emissions.

For example, the sum of the hearing system rating of 4 (M4) and a telephone rating of 3 (M3) would result in a combined rating of 7. Any combined rating that equals at least 5 would provide "normal use"; a combined rating of 6 or greater would indicate "excellent performance." The immunity of this hearing system is at least M4. The equipment performance measurements, categories and system classifications are based upon the best information available but cannot guarantee that all users will be satisfied.

The performance of the individual hearing aids may vary with individual mobile phones. Therefore, please try this hearing aid with your mobile phone or, if you are purchasing a new phone, please be sure to try it with your hearing aid prior to purchase. For additional guidance, please ask your hearing healthcare professional for the booklet entitled.

## Technical Information.

Changes or modifications not expressly approved by the party responsible for compliance could void the users authority to operate the equipment.

### Radio information of your wireless hearing

Antenna type	Resonant loop antenna
Operation frequency	2.402–2.480 GHz
Modulation	GFSK, Pi/4 DPSK
Radiated power	<1mW (EIRP)

### Bluetooth®

Range	~ 1m (class 3)
Bluetooth	Version 4.2
Profiles supported	HFP (Hands-free profile)

D **Operating Conditions:**

This device is designed such that it functions without problems or restrictions if used as intended, unless otherwise noted in these user guides.

F **Transport and Storage Conditions:**

Temperature:  $-20^{\circ}$  to  $+60^{\circ}$  Celsius ( $-4^{\circ}$  to  $+140^{\circ}$  Fahrenheit).

E Humidity: Up to 90% (non condensing).

### Disposal Information.

Within the European Union, the marked equipment is covered by „Directive 2002/96/EC of the European Parliament and of the Council of 27 January 2003 on waste electrical and electronical equipment“. Amended by „Directive 2003/108/EC“ (WEEE).

Recycle the device and batteries according to national regulations.

### Compliance Information.

#### Declaration of Conformity

Hansaton hereby declares that this Hansaton product meets the requirements of the Medical Device Directive 93/42/EEC as well as the Radio Equipment Directive 2014/53/EU.

The full text of the Declaration of Conformity can be obtained from the manufacturer.

#### Australia/New Zealand:

**R-NZ**

Indicates a device's compliance with applicable Radio Spectrum Management's (RSM) and Australian Communications and Media Authority (ACMA) regulatory arrangements for the legal sale in New Zealand and Australia. The compliance label R-NZ is for radio products supplied in the New Zealand market under conformity level A1.

#### Notice 1

This device complies with Part 15 of the FCC Rules and with RSS- 210 of Industry Canada. Operation is subject to the following two conditions:

- 1) this device may not cause harmful interference, and
- 2) this device must accept any interference received, including interference that may cause undesired operation.

#### Notice 2

Changes or modifications made to this device not expressly approved by Hansaton may void the FCC authorization to operate this device.

**Notice 3**

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and ICES-003 of Industry Canada. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and receiver
- Connect the device into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help.

**FDA Information (Valid for USA only)****Warning to hearing aid dispensers (to comply with the U.S. Food and Drug Administration regulations)**

A hearing aid dispenser should advise a prospective hearing aid user to consult promptly with a licensed physician (preferably an ear specialist) before dispensing a hearing aid if the hearing aid dispenser determines through inquiry, actual observation, or review of any other available information concerning the prospective user, that the prospective user has any of the following conditions:

- (i) Visible congenital or traumatic deformity of the ear.
- (ii) History of active drainage from the ear within the previous 90 days.
- (iii) History of sudden or rapidly progressive hearing loss within the previous 90 days.
- (iv) Acute or chronic dizziness.
- (v) Unilateral hearing loss of sudden or recent onset within the previous 90 days.
- (vi) Audiometric air-bone gap equal to or greater than 15 decibels at 500 hertz (Hz), 1,000 Hz, and 2,000 Hz.
- (vii) Visible evidence of significant cerumen accumulation or a foreign body in the ear canal.
- (viii) Pain or discomfort in the ear. Special care should be exercised in selecting and fitting a hearing aid whose maximum sound pressure level exceeds 132 decibels because there may be risk of impairing the remaining hearing of the hearing aid user. [This provision is required only for those hearing aids with a maximum sound pressure capability greater than 132 decibels (dB).]

## D **Important Notice for Prospective Hearing Aid Users.**

EN Good health practice requires that a person with a hearing loss have a medical evaluation by a licensed physician (preferably a physician who specializes in diseases of the ear) before purchasing a hearing aid.

F Licensed physicians who specialize in diseases of the ear are often referred to as otolaryngologists, otologists or otorhinolaryngologists. The purpose of medical evaluation is to assure that all medically treatable conditions that may affect hearing are identified and treated before the hearing aid is purchased. Following the medical evaluation, the physician will give you a written statement that states that your hearing loss has been medically evaluated and that you may be considered a candidate for a hearing aid. The physician will refer you to an audiologist or a hearing aid dispenser, as appropriate, for a hearing aid evaluation. The audiologist or hearing aid dispenser will conduct a hearing aid evaluation to assess your ability to hear with and without a hearing aid. The hearing aid evaluation will enable the audiologist or dispenser to select and fit a hearing aid to your individual needs. If you have reservations about your ability to adapt to amplification, you should inquire about the availability of a trial-rental or purchase-option program. Many hearing aid dispensers now offer programs that permit you to wear a hearing aid for a period of time for a nominal fee after which you may decide if you want to purchase the hearing aid. Federal law restricts the sale of hearing aids to those individuals who have obtained a medical evaluation from a licensed physician. Federal law permits a fully informed adult to sign a waiver statement declining the medical evaluation for religious or personal beliefs that preclude consultation with a physician. The exercise of such a waiver is not in your best health interest and its use is strongly discouraged.

### **Children With Hearing Loss.**

In addition to seeing a physician for a medical evaluation, a child with a hearing loss should be directed to an audiologist for evaluation and rehabilitation since hearing loss may cause problems in language development and the educational and social growth of a child. An audiologist is qualified by training and experience to assist in the evaluation and rehabilitation of a child with a hearing loss.

### **Cell phone.**

Some hearing aid users have reported a buzzing sound in their hearing aids when they are using cell phones – indicating that the cell phone and hearing aid may not be compatible. According to the ANSI C63.19 standard (ANSI C63.19-2011 American National Standard Methods of Measurement of Compatibility Between Wireless Communications Devices and Hearing Aids), the compatibility of a particular hearing aid

and cell phone can be predicted by adding the rating for the hearing aid immunity to the rating for the cell phone emissions. For example, the sum of the hearing aid rating of 2 (M2) and a telephone rating of 3 (M3) would result in a combined rating of 5. Any combined rating that equals at least 5 would provide "normal use": a combined rating of 6 or greater would indicate "excellent performance". The immunity of this hearing aid is at least M2/T2. The equipment performance, measurements, categories and system classifications are based upon the best information available but cannot guarantee that all users will be satisfied.

D

GB

F

E

Note: The performance of the individual hearing aids may vary with individual cell phones. Therefore, please try this hearing aid with your cell phone or, if you are purchasing a new phone, please be sure to try it with your hearing aid prior to purchase. For additional guidance, please ask your hearing healthcare professional.

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Hansaton Akustik GmbH